

SINCE 1978



OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED
An ISO 10002:2004 CERTIFIED

CE Certified

Collapsible tube filling and sealing equipment



Stroke type, barrel and piston arrangement with stroke setting for various capacity of filling of tubes (from 5 gms to 50 gms.) detachable nozzles provided for each separate capacity. All contact parts are of stainless steel also stainless steel holder. Built on cast iron body, hand lever operated.

Collapsible tube filling and sealing equipment