

SINCE 1978



# OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED  
An ISO 10002:2004 CERTIFIED

**CE Certified**

## Cooling Water Bath



Cooling Water Bath

Double walled, Inner Chamber S.S. Outer M.S. Duly Powder Coated / S.S.  
PUF Insulated  
Temp. Range 5°C to 60°C.  
Cooling Compressor.  
Copper Cooling Coils Welded throughout on the Outside of the S.S. Chamber.  
Body Mounted on sturdy Angle Iron Frame with Caster wheels for easy mobility.  
Temp. Is controlled by Digital Temp. Controller.  
Provided with circulating Pump.

Working Size : 12" x 12" x 12"

### Optional :

Microprocessor Based P.I.D. controller