

SINCE 1978



OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED
An ISO 10002:2004 CERTIFIED

CE Certified

De-ionizer (de-mineralizer)



De-ionizer (de-mineralizer)

It consist of an inbuilt prefilter & very high exchanging capacity non-corrosive CATION & ANION Resin, columns prefilled with imported Resins for efficient results having treated water parameters. Fitted with a sturdy PVC multi control valve & can be connected direct to water taps. Fitted with analog / digital conductivity meter. Supplied with one plastic chemical proof regeneration tank & the complete unit is housed in a painted M.S. trolley for easy mobility.

FEATURES:

- Non corrosive Regeneration Tank
- In Built Pre filter
- High capacity Cation and Anion exchange resins.
- Upflow treatment.
- Easily portable - metal stand and trolley
- Can be connected direct to Water Tap.
- Simple Periodic regeneration - negligible operating costs.
- Digital model is fitted with digital conductivity meter.
- Treated water will have the conductivity less than 10 ms/cm and pH of 7.5 - 9
- Standard models available in 50, 70, 100, 200, 300 Ltrs / Hr.

Also available in FOUR bed and mixed bed models

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: info@opticstechnology.in, opticstechnology@gmail.com

Website: www.opticstechnology.in