

SINCE 1978



OPTICS TECHNOLOGY



Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED
An ISO 10002:2004 CERTIFIED

CE Certified

Digital Colorimeter



Digital Colorimeter

- Stable, direct readout of absorbance
- Compact & easy to operate
- Rugged, sturdy & reliable
- Thumb wheel selection of filters
- Pre - focused light source
- long life photo diode detector
- 1 ml sample size

WAVE LENGTH RANGE	400 to 700 nm
RESOLUTION	0.01
ABSORBANCE	0 – 1.99A
PHOTO DETECTOR	Silicon Photodiode / Photo Cell
DISPLAY	2½ digit LED
LIGHT SOURCE	6.8V, 0.3Amp. Tungsten lamp
SAMPLE SYSTEM	10mm path length matched glass test tubes
FILTERS	8 Filters : 400, 420, 480, 500, 520, 540, 620, 680nm
POWER SUPPLY	230V AC ±10%, 50Hz

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: info@opticstechnology.in, opticstechnology@gmail.com

Website: www.opticstechnology.in