

SINCE 1978



# OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED  
An ISO 10002:2004 CERTIFIED

**CE Certified**

## Digital Potentiometer



Digital Potentiometer

For Highly Accurate & Reliable Potentiometer Measurements in the Laboratory

- Highly Stable & Accurate
- 4 Electrodes For Different Titrations
- 3 ½ Digit LED Display
- Rugged & Economical

RANGE	0 to ± 199.9mV / 0 to ± 1999 mV
RESOLUTION	0.1mV
REPEATABILITY	±1mV
ACCURACY	±1mV ±1 digit
INPUT IMPEDANCE	> 1012 ohms
OPERATING TEMPERATURE	10° C to 45° C
DISPLAY	3 ½ Digit 0.5" seven segment LED display with auto polarity indication

### STANDARD ACCESSORIES

Platinum Electrode	1 No.
Reference Electrode	1 No.
Glass Electrode	1 No.
Silver Electrode	1 No.
Electrode Stand & Clamp	1 No.

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: [info@opticstechnology.in](mailto:info@opticstechnology.in), [opticstechnology@gmail.com](mailto:opticstechnology@gmail.com)

Website: [www.opticstechnology.in](http://www.opticstechnology.in)