

SINCE 1978



# OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED  
An ISO 10002:2004 CERTIFIED

**CE Certified**

## Digital Whiteness Meter



Digital Whiteness Meter

Whiteness Meter based on the principle of photoelectric effect & D65 light source (R457) as a illumination condition. Meanwhile except solution or liquid samples, this instrument can be widely used for measuring the whiteness value of any substances Such as: powder, board, paper, paste, sheet, clay etc. This instrument is used widely f or measuring the whiteness value of any substances such as powder, board, paper, plastic, sheet, toothpaste, fabric etc.

Measurement Range	0 - 100%
Accuracy	0.1
Repeatability	0.2
Power supply	220V $\pm$ 10%

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: [info@opticstechnology.in](mailto:info@opticstechnology.in), [opticstechnology@gmail.com](mailto:opticstechnology@gmail.com)

Website: [www.opticstechnology.in](http://www.opticstechnology.in)