

SINCE 1978



# OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED  
An ISO 10002:2004 CERTIFIED

**CE Certified**

## Floculator / Jar Test Apparatus



Floculator / Jar Test Apparatus

Jar Test Apparatus allows efficient and economical flocculation matter and collides in general and also has been specially designed for use in water treatments to correctly estimate the dosing of alum and such other coagulants for treatment of Water sewage and floc formation test in Water Treatment plants.

Fitted with D.C. motor but to work on 220 V AC. Variable speed from 20 to 100 RPM approx. Alongwith illumination

Capacity of flasks 1 Ltr / 2 Ltrs

For 2 Stirrers  
For 4 Stirrers  
For 6 Stirrers  
For 8 stirrers

### OPTIONAL

Mechanical Timer

**OR**

Digital Timer

With Digital RPM Indicator (Extra)

Attachment for feeding dose to all the stirring units simultaneously

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: [info@opticstechnology.in](mailto:info@opticstechnology.in), [opticstechnology@gmail.com](mailto:opticstechnology@gmail.com)

Website: [www.opticstechnology.in](http://www.opticstechnology.in)