

SINCE 1978



OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED
An ISO 10002:2004 CERTIFIED

CE Certified

Fumigator



Fumigator

Fumigators for Pharmaceutical, Food & Dairy, Biotech industries generate high density fog with droplets size of 5-15 microns, which helps in providing maximum thrust, deeper penetration & wide spread of aerosols. Some of the common application areas it is used in this industry include injectable, vaccines, tablets, parenterals, IV fluids, medical devices, bulk drug manufacturing areas, QC & QA areas, microbiology lab, R&D center, clinical drug trial centre, tissue culture laboratories, pathology laboratories and many more.

OPTIONAL

(a) Humidistat (To Control Humidity)

OR

Digital Humidity Controller

(b) Mechanical Timer (0-60 min)

OR

Digital Timer

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: info@opticstechnology.in, opticstechnology@gmail.com

Website: www.opticstechnology.in