

SINCE 1978



OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED
An ISO 10002:2004 CERTIFIED

CE Certified

Gerber Centrifuge



Gerber Centrifuge



Gerber Centrifuge

Electrically Operated Gerber Centrifuge Machine fitted with Heavy duty Electric Motor housed in heavy cast iron body, Timer for 0-5 Min. & with Clamps, dynamically balanced body.

Models available:-

- (8 Test) With Brake & Timer
- (12 Test) With Brake & Timer
- (24 Test) With Brake & Timer

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: info@opticstechnology.in, opticstechnology@gmail.com

Website: www.opticstechnology.in