

SINCE 1978



# OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED  
An ISO 10002:2004 CERTIFIED

**CE Certified**

## Germinator Seed



Germinator Seed

Most suitable for conducting various experiments on variety of seeds under different conditions of temperature and humidity. Atmospheric conditions are created within the chamber.

- Temperature from 5oC to 50o C  $\pm$  1oC Humidity range Ambient to 95%  $\pm$  3%.
- Temp. & Humidity is controlled by Digital Temp. cum Humidity Controller.
- Each chamber has 14 trays.
- The Instrument has stainless steel water Reservoir at back for humidity.
- Double walled inside made of Stainless Steel and Outside made of Mild Steel with Powder Coating.
- External lighting arrangement by tube lights and Mechanical Cyclic Timer 0-24 hrs.
- Work to on 220/230V A.C.

Dimension of each chamber (Inside)  
(WXDXH) 600 X 600 X 900 mm (Approx.)  
24" X 24" X 36" (Approx.)

### OPTIONAL:-

- A) Automatic cyclic timer (0-24 Hrs.)
- B) Dual setting digital temperature controller for setting two different temperatures for day and night.
- C) Automatic voltage stabilizer
- D) Microprocessor based Temp. & Humi. Controller



**NOTE: DUAL CHAMBER MODEL & OUTER STAINLESS STEEL MODELS ARE ALSO AVAILABLE TOUCH SCREEN DISPLAY MODELS WITH DATA LOGGING OPTIONA ARE ALSO AVAILABLE ON REQUEST.**

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: [info@optics technology.in](mailto:info@optics technology.in), [optics technology@gmail.com](mailto:optics technology@gmail.com)

Website: [www.opticstechnology.in](http://www.opticstechnology.in)