

SINCE 1978



OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED
An ISO 10002:2004 CERTIFIED

CE Certified

Gyratory Sieve Shaker



Gyratory Sieve Shaker

FOR 8", 12" & 18" TEST SIEVES Gyratory Sieve Shaker is driven by ¼ H.P. Electric motor through a reduction gear, can carry upto SIX, 6" OR 8" /12" /18" Test Sieves with one set of LID AND Pan . In this Gyratory motion of the sieve shaker there is an upward and downward movement but is inclined from the vertical axis and the direction or inclination changes progressively in clockwise direction having frequency of about 270 vibration per minutes. Provided with or without timer as per requirement.

SIEVES IN DIFFERENT MESH SIZES OF DIA 8", STAINLESS STEEL WIRE MESH IN BRASS / STAINLESS STEEL FRAME WITH OR WITHOUT CALIBRATION ARE ALSO AVAILABLE ON REQUEST.



2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: info@opticstechnology.in, opticstechnology@gmail.com

Website: www.opticstechnology.in