

SINCE 1978



OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED
An ISO 10002:2004 CERTIFIED

CE Certified

Heavy Duty Willey Mill



Heavy Duty Willey Mill

SPECIFICATIONS:-

Effective and Continuous crushing, leading to increase production and faster return on investments.

Heavy Duly Willey Mill for Continuous working.

Suitable for Wet-Pulpy and Dry crushing.

Compact Size. (Length 56 cm, Width 38 cm, Height 69 cm) No wastage of product and no pollution hazards. The Willey Mill grinds almost all Spices, green chillies, fresh ginger, garlic, salt, sugar, coffee, fresh fruits, fruit pulp, vegetables, tobacco, herbs, beetle nuts, herbs roots, leaves, bark of trees, oil cake, cattle and poultry feeds, icing sugar, costarred powder, instant food, food colors and allied food production.

Provided with Set of 5 Nos. perforated round hole screens, 1 No. Brush for cleaning of machine

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: info@opticstechnology.in, opticstechnology@gmail.com

Website: www.opticstechnology.in