

SINCE 1978



# OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED  
An ISO 10002:2004 CERTIFIED

**CE Certified**

## Horizontal Gel Electrophoresis



### DESCRIPTION

Electrophoresis of nucleic acid (dna/rna) in agarose gel has become an essential step in every walk of biotechnology today. High buffer-gel ratio is kept in all our systems whereby current flow is minimized; temperature rise and pH, difference are also kept at minimum levels. This results in fast and good resolution of bands.

### FEATURES:

- Basic unit with platinum electrode assembly.
- Comb 2mm thick
- 10 well and thick comb stand
- Lid
- Connecting cord

Homogenizer Emulsifier

### DIGITAL ELECTROPHORESIS POWER SUPPLY

It delivers stabilized output 0-500 Volts DC and 0-100mA. Individually and continuously adjustment with three figure digital display of volts milliamps irrespective of the mains fluctuation and variation in paper resistance. It is completely solid state to deliver constant voltage and current as desired.



DIGITAL ELECTROPHORESIS POWER SUPPLY

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: [info@opticstechnology.in](mailto:info@opticstechnology.in), [opticstechnology@gmail.com](mailto:opticstechnology@gmail.com)

Website: [www.opticstechnology.in](http://www.opticstechnology.in)