

SINCE 1978



OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED
An ISO 10002:2004 CERTIFIED

CE Certified

LED based Slim X-Ray Viewer



LED based Slim X-Ray Viewer

Ultra slim, compact, solid with Backlight LED's based panel

- Long Life 50,000 Hrs. (Approx.)
- Uniform Light all over the surface area
- No Heat emission
- Wall mounted
- Power consumption very low 12V (DC), 2 A
- Screen size 300 x 350 mm & Thickness 19mm

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: info@opticstechnology.in, opticstechnology@gmail.com

Website: www.opticstechnology.in