

SINCE 1978



OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED
An ISO 10002:2004 CERTIFIED

CE Certified

Microprocessor Based pH Meter



Microprocessor Based pH Meter

Display	16 x 2 Alphanumeric LCD display	
Range	pH: 0-14 pH,mV: 0 to ±1999	Temp:0 to 100°C
Accuracy	pH: ± 0.01pH,mV: ±1mV	0.1°C
Calibration	2 Point OR 3 point	
Temp. Compensation	Auto: 0 to 100°C,	Manual: 0 to 100°C
Slope Control	80 to 120%	
Power	220 / 230V AC, 50 Hz	
Sensor	pH electrode	
Auto Buffer recognition, Storage for 90 samples, Printer Attachment Facility (Optional)		

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: info@opticstechnology.in, opticstechnology@gmail.com

Website: www.opticstechnology.in