

SINCE 1978



OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED
An ISO 10002:2004 CERTIFIED

CE Certified

Oil Bath



Oil Bath

- Double Wall Construction
 - Inner Chamber SS-304 & Outer M.S. Duly Powder coated
 - with Heavy Duty Argon Welded tank
 - Temperature Range Amb. +5 to 200°C
 - Thermostatic control of temperature
 - Stirrer for uniform temperature inside the working chamber
- Size 250 X 250 X 250 mm

OPTIONAL

Digital Temperature Controller-cum-Indicator
OR
Microprocessor Based PID Controller

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: info@opticstechnology.in, opticstechnology@gmail.com

Website: www.opticstechnology.in