

SINCE 1978



# OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED  
An ISO 10002:2004 CERTIFIED

**CE Certified**

## Portable Autoclave (Pressure Cooker Type)



Electrically heated, complete with stand single walled with Stainless Steel Body along with essential parts like pressure gauge, safety valve etc. Outer, Inner & Lid is made of stainless steel. The chamber is absolutely leak proof. It is provided with safety valve, pressure gauge, steam release valve, and vacuum release valve.

**Available sizes:**

300 x 300 mm

Cap. 20 Ltrs. Approx.

Portable Autoclave (Pressure Cooker Type)

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: [info@opticstechnology.in](mailto:info@opticstechnology.in), [opticstechnology@gmail.com](mailto:opticstechnology@gmail.com)

Website: [www.opticstechnology.in](http://www.opticstechnology.in)