

SINCE 1978



# OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED  
An ISO 10002:2004 CERTIFIED

**CE Certified**

## RECEIPROCATED KAHN SHAKER



RECEIPROCATED KAHN SHAKER

D.C. Motor special system shaker

Works on 220 Volts A.C.

For variable and smooth function speed control by speed regulator.

It is general purpose machine and found useful for mixing or agitation ingredients, solution in flasks, bottles, beakers and kahn racks etc.

Upper tray covered with rubber floor and sides of the tray are slotted to take parallel and cross bars cover with rubber. The bars are tightened with washers and fly nuts which can be adjusted as per the requirement.

### PLATFORM SIZE :

(A) 20" x 13" x 4"

(B) 31" x 17" x 4"

(C) 36" x 18" x 4"

### Optional :

(A) Mechanical Timer

**OR**

(B) Digital Timer

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: [info@opticstechnology.in](mailto:info@opticstechnology.in), [opticstechnology@gmail.com](mailto:opticstechnology@gmail.com)

Website: [www.opticstechnology.in](http://www.opticstechnology.in)