

SINCE 1978



OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED
An ISO 10002:2004 CERTIFIED

CE Certified

Wrist Action Shaker



Wrist Action Shaker

The shaking machine fitted with D.C. MOTOR to work on 220 VOLT A/C 50 CYCLES to work for longer period smoothly. Flasks mounted spindle to oscillate through a small amplitude sturdy eccentric mechanism provides effective mechanical shaking . The shaking machine constructed to hold 100 ml to 1000 ml. Conical flasks. The shaking speed action is controlled by ELECTRONIC regulator . 100 R.P.M TO 800 R.P.M approx

HOLDING CAPACITY

- (A) To hold 4 Flasks Cap. 100 ml. to 1000 ml. by side clamps
(B) To hold 8 Flasks Cap. 100 ml. to 1000 ml. by side clamps

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: info@opticstechnology.in, opticstechnology@gmail.com

Website: www.opticstechnology.in