

SINCE 1978



# OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED  
An ISO 10002:2004 CERTIFIED

**CE Certified**

## Slide Warmer (For I.V.F. Labs)



Slide Warmer (For I.V.F. Labs)

Digital units include built-in temperature sensing probe Compact design saves valuable bench space.  
Powder-coated steel body construction ensures durability.  
Digital units utilize a Digital Temperature Controller and easy-to-read LED display for accurate and reproducible temperature control.  
Temperature Range: ambient +5°C to 50°C.  
Petri plate warmer, S.S. body, Anodized Aluminium Heating Plate,  
Temperature controlling by Digital Temperature controller, Size 300 x 100 x 40 mm Approx

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: [info@opticstechnology.in](mailto:info@opticstechnology.in), [opticstechnology@gmail.com](mailto:opticstechnology@gmail.com)

Website: [www.opticstechnology.in](http://www.opticstechnology.in)