

SINCE 1978



OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED
An ISO 10002:2004 CERTIFIED

CE Certified

Soxhlet Extraction Apparatus



Soxhlet Extraction Apparatus

The complete soxhlet apparatus consists of solvent reservoir extraction assembly with 6 extraction places with concentric metal rings of various diameter to suit 100ml flasks to 250ml flasks and built in steam generator.

Supplied with stand, clamp holder soxhlet extraction (condenser) glass and extraction flask of 60 ml volume.

Mantle type heaters are provided with individual energy regulators.

Available for 3 unit test OR 6 unit Test.

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: info@opticstechnology.in, opticstechnology@gmail.com

Website: www.opticstechnology.in