

SINCE 1978



OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED
An ISO 10002:2004 CERTIFIED

CE Certified

Table Top Sieve Shaker



Table Top Sieve Shaker

(FOR 8" & 12" TEST SIEVES) (Soundless , Vibration type) Fitted with F.H.P Electric Motor. It can accommodate 6 Test sieves of 8" Dia with Lid and pan. Speed can be regulated through variable speed regulator to work on 220 volts A.C (without Test Sieves) speed regulator to work on 220 volts A.C

SIEVES IN DIFFERENT MESH SIZES OF DIA 8", STAINLESS STEEL WIRE MESH IN BRASS / STAINLESS STEEL FRAME WITH OR WITHOUT CALIBRATION ARE ALSO AVAILABLE ON REQUEST.

SIEVES IN DIFFERENT MESH SIZES OF DIA 8", STAINLESS STEEL WIRE MESH IN BRASS / STAINLESS STEEL FRAME WITH OR WITHOUT CALIBRATION ARE ALSO AVAILABLE ON REQUEST.



2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: info@opticstechnology.in, opticstechnology@gmail.com

Website: www.opticstechnology.in