

SINCE 1978



OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED
An ISO 10002:2004 CERTIFIED

CE Certified

Tablet Friability Test Apparatus



Tablet Friability Test Apparatus

The apparatus is designed to determine the abrasion rate and impact hardness of tablets. One/Two transparent acrylic drums rotated by superior quality constant speed sturdy motor, mounted in painted sheet metal cabinet with chemical resistant epoxy coating.

SPECIAL FEATURES: Realistic abrasion rate is obtained by pre-set revolutions through a pilot lamp as compared to time setting, where voltage fluctuation may change number of revolution in pre-set time.

SPECIFICATIONS:

Drum Diameter :285mm
Tablet Fall Distance:140mm
Drive :1/8 HP
Power Phase.:AC/230 volt Single

Available in SINGLE AND DOUBLE DRUM MODELS WITH OUTER BODY OF CONSTRUCTION BEING M.S.DULY POWDER COATED BODY OR S.S. 304.

Optional :

Digital Revolution Counter.
Spare Drum.



2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: info@opticstechnology.in, opticstechnology@gmail.com

Website: www.opticstechnology.in