

SINCE 1978



OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED
An ISO 10002:2004 CERTIFIED

CE Certified

THIN LAYER CHROMATOGRAPHY (T.L.C. APPARATUS)



Thin Layer applicator is specially constructed to give dependable service in drawing thin uniform layer variable form 0 to 2 mm. (Minimum thickness attainable 0.1 mm.)
Accurate design and skilled workmanship of the applicator gives no trouble in coating number of plates at a time. It is made of Aluminum finished in proper shape.

This kit consists of the following:

1. Applicator
2. Drawing Board
3. Drying Rack
4. Spotting Guide
5. Developing Chamber
6. Glass Pipette
7. Silica Gel
8. Glass Plates
9. Glass Jar (Optional)

THIN LAYER CHROMATOGRAPHY
(T.L.C. APPARATUS)

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: info@opticstechnology.in, opticstechnology@gmail.com

Website: www.opticstechnology.in