

SINCE 1978



OPTICS TECHNOLOGY

Manufacturer & Exporters of : **LABORATORY TESTING INSTRUMENTS**



An ISO 9001:2008 CERTIFIED
An ISO 10002:2004 CERTIFIED

CE Certified

U.V. Transilluminator



U.V. Transilluminator

Transilluminators find wide applications for DNA analysis. All models have high output UV tubes, which have an average life expectancy of 5000 hours. The UV filters are polished to give blemish free viewing of the gel.

COLLEGE MODEL TRANSILLUMINATOR WITH FILTER

Wavelength	:	312nm
Fixed No.of	:	Tubes 6
Wattage	:	8 Watts
Detection limit	:	200 ng
Filter area	:	20X20 cm
Body	:	MS with powder coated

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: info@opticstechnology.in, opticstechnology@gmail.com

Website: www.opticstechnology.in