

SINCE 1978



OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED
An ISO 10002:2004 CERTIFIED

CE Certified

Water Bath (Cooling)



Water Bath (Cooling)

Double walled, Inner Chamber S.S. Outer M.S. Duly Powder Coated / S.S.
PUF Insulated
Temp. Range 5°C to 60°C.
Cooling Compressor.
Copper Cooling Coils Welded throughout on the Outside of the S.S. Chamber.
Body Mounted on sturdy Angle Iron Frame with Caster wheels for easy mobility.
Temp. Is controlled by Digital Temp. Controller.
Provided with circulating Pump.

Working Size : 12" x 12" x 12"

Optional :

Microprocessor Based P.I.D. controller

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: info@opticstechnology.in, opticstechnology@gmail.com

Website: www.opticstechnology.in