

SINCE 1978



OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED
An ISO 10002:2004 CERTIFIED

CE Certified

Water Bath Incubator Shaker



(Metabolic Shaking Incubator)

Temperature range from room temperature to 100°C.
Temp. is Controlled by Thermostat / Digital Temp. Controller / PID with an accuracy of $\pm 0.5^{\circ}\text{C}$.
Double walled inside made of stainless steel 304 Quality and outside Mild Steel sheet painted in epoxy powder coating/S.S.
Oscillating tray is riding on ball rollers.
Shaking speed range 40 to 140 rpm.
To work on 220/230 volts A.C. Supply.

SHAKING TRAY SIZE: INSIDE
TRAY MADE OF S.S. 304 QUALITY)
Models :

Water Bath Incubator Shaker

A)	275x275x75 mm.	6 Ltrs.
B)	275x275x150mm.	12 Ltrs.
C)	405x300x150mm	18 Ltrs.
Cooling Models With Temp. Range 5 - 50°C are also available		



2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)
Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890
E-mail: info@opticstechnology.in, opticstechnology@gmail.com
Website: www.opticstechnology.in