

SINCE 1978



OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED
An ISO 10002:2004 CERTIFIED

CE Certified

Water Distillation (S.S.)



Water Distillation (S.S.)

WALL MOUNTED - Heavy Duty Argon Welded

- The unit is recommended for making pyrogen free distilled water for routine laboratory requirements.
- Made out of Stainless Steel (SS-304 grade) sheet, Heavy Duty with Argon Welding.
- A wall- mounting clamp is provided with the unit.
- S.S.-304 grade Argon Arc Welded joints.
- Supplied with cord and plug to work on 220/230 volts A.C. supply.
- Heaters are of immersion type
- The Still is provided with bracket (to be fixed on wall) on which the still can rest.

Capacity in Ltrs.

- 2 Litre/hrs
- 4 Litre/hrs
- 6 Litre/hrs
- 8 Litre/hrs
- 10 Litre/hrs

Optional:

Water Level Cut off System for safety of heaters in case of No / Insufficient water.



TABLE PATTERN -

- Double Walled Construction.
- Inner & Outer S.S. 304. The gap in between being filled with glass wool.
- The condenser is made of S.S. to work on 220/230 volts
- Delux Model Heavy Duty with Argon Welding

Capacity in Ltrs.

- 2 Litre/hrs
- 4 Litre/hrs
- 6 Litre/hrs
- 8 Litre/hrs
- 10 Litre/hrs

Optional:

Water Level Cut off System for safety of heaters in case of No / Insufficient water.

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: info@opticstechnology.in, opticstechnology@gmail.com

Website: www.opticstechnology.in