

SINCE 1978



# OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED  
An ISO 10002:2004 CERTIFIED

**CE Certified**

## Willey Mill



Willey Mill

INTERMEDIATE MODEL "ARTHUR THOMAS" TYPE Electrically operated to work on 220 volts A.C Recommended for grinding of dry plant samples. Also for a great variety of commercial dry materials. Unit having two stationary steel blades and rotary with steel edges provided in the cast iron chamber.

Willey mill is fitted on substantial cast iron base.

### WORKING CHAMBER

Size 40 X 25 mm  
Size 65 X 25 mm  
Size 75 x 25 mm  
Size 100 x 50 mm  
Size 200 x 75 mm

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: [info@opticstechnology.in](mailto:info@opticstechnology.in), [opticstechnology@gmail.com](mailto:opticstechnology@gmail.com)

Website: [www.opticstechnology.in](http://www.opticstechnology.in)