



OPTICS TECHNOLOGY

SINCE 1978

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2018 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

ABEL'S FLASH POINT APPARATUS

CE



OT-ASP-248

Electrically Heated model with voltage Variac for control. Suitable for determination of closed cup flash point determination of Petroleum and petroleum products having flash points below 70° C.

Optional :

- Digital temperature indicator
- Electrically Driver Stirrer with a FHP Motor having flexible shaft.
- Oil Cup Thermometer (-35°C to 70°C).
- Water Bath Thermometer (-30°C to 80°C)

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in

opticstechnology@gmail.com

Website : www.opticstechnology.in