



OPTICS TECHNOLOGY

SINCE 1978

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2018 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

COD DIGESTOR CE



OT-COD-145

Heating block type Temperature range from 5°C above ambient to 150°C \pm 1°C.

Designed with glass reactions vessels and air condensers.

Supplied with stand to accommodate vessels Operate on 230 Volts, 50 Hz, single phase.

With Digital Temperature controller cum digital timer along with glass parts.

Model No. **No. of Test**

OT-COD-145/A 06 Sample

OT-COD-145/B 15 Sample

OT-COD-145/C 20 Sample

OT-COD-145/D 25 Sample

Or any other Size on Request.

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in

opticstechnology@gmail.com

Website : www.opticstechnology.in