



OPTICS TECHNOLOGY



SINCE 1978

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2018 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

DIGITAL REFRACTOMETER (PORTABLE)

(Digital Indication of Brix & Refractive Index CRI)



OT-DR-05

1. Waterproof
2. Automatic temperature compensation (ATC)
3. Temperature display
4. High/Low range limit alerts
5. Auto shut-off and a low battery indicator
6. Built - in Brix % and nD scales
7. Powered by Battery or AC

Technical Details:

Test Item	Measuring Range	Resolution	Accuracy
Brix	0.0-95%	0.1	±0.1%
RI(nD)	1.3330-1.5318	0.0001	±0.0002
Temp	5-80°C	0.1°C	±1°C

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in,

opticstechnology@gmail.com,

Website : www.opticstechnology.in