



# OPTICS TECHNOLOGY



SINCE 1978

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010  
An ISO 10002:2018 CERTIFIED BS EN 12469:2000  
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022


## HAND REFRACTOMETER

### (POCKET REFRACTOMETER)



OT-HR-01

We are sole manufacturer of Refractometer with "BRASS BODY" in India



Model No.

OT-HR-01/32

OT-HR-01/62

OT-HR-01/92

OT-HR-01/55

OT-HR-01/82

Range :

(A) (0 - 32%)

(B) (28 - 62%)

(C) (58 - 92%)

(D) (0 - 55%)

(E) (45 - 82%)

OT-HR-01/80

(F) (0 - 90%)

New Specially Designed Full Range **0-90%**

We are sole manufacturer of Refractometer with "BRASS BODY" in India

### (Brass Body)

Hand Refractometer is a light-based instrument for field and laboratory uses for estimating sugar contents, Brix in T.S.S.

"Optics Technology" make Hand Refractometer (Brix Meter) is specially designed for measuring sugar content in an aqueous solution, soft drinks, sugar solutions, cutting liquid, fruit juice or milk and beverages.

Easy to use Cushioned with soft & comfortable non-slip rubber grip, Adjustable focus and clearly defined scale, with easy-to-read measurements.

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: [info@opticstechnology.in](mailto:info@opticstechnology.in),

[opticstechnology@gmail.com](mailto:opticstechnology@gmail.com),

Website : [www.opticstechnology.in](http://www.opticstechnology.in)