



OPTICS TECHNOLOGY



SINCE 1978

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2018 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

MICROPROCESSOR BASED pH METER (Table Top) CE



OT-pH-105-M

Display	: 16 x 2 Alphanumeric LCD display
Range	: pH: 0-14 pH, mV: 0 to ±1999, Temp: 0 to 100°C
Accuracy	: pH: ± 0.001pH, mV: ±1mV: 0.1°C
Calibration	: 2 Point OR 3 Point OR 5 Point
Temp. Compensation	: Auto: 0 to 100°C, Manual: 0 to 100°C
Slope Control	: 80 to 120%
Power	: 220 / 230V AC, 50 Hz
Sensor	: pH electrode (Glass)

Auto Buffer recognition, Storage for 90 samples, Printer Attachment Facility (Optional)

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in

opticstechnology@gmail.com

Website : www.opticstechnology.in