



SINCE 1978

OPTICS TECHNOLOGY



Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2018 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

RESEARCH BINOCULAR MICROSCOPE

CE



OT-BMS-161 R

- Binocular Microscope standard set complete with 6V - 20 W halogen light illuminator
- With anti-fungus Achromatic Objectives 4x, 10x, 40x(spring) & 100x (oil, immersion)
- Paired eyepieces wide field 10x and illumination through SMPS circuit
- Abbe condenser 1.25 with aperture iris diaphragm
- Coaxial focus, tension adjustment on right.
- Quadruple revolving nose piece.

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in

opticstechnology@gmail.com

Website : www.opticstechnology.in