



SINCE 1978

# OPTICS TECHNOLOGY

**Manufacturer & Exporters of : LAB TESTING INSTRUMENTS**



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010  
An ISO 10002:2018 CERTIFIED BS EN 12469:2000  
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

**VACUUM OVEN**

**ROUND**

CE



**OT-VO-136 MS**



**OT-VO-136 SS**

Inside chamber made of S.S. 304 Outside Mild Steel / S.S. 304 with Toughened Glass Door. It is provided with 1 or 2 No. of S.S. Trays.

Temp. Range Ambient Temp. +5°C to 130°C. Temp. Controlled by microprocessor based PID controller. Vacuum Pressure indication by vacuum gauge.

| Model No.         | Dia mm x Depth mm           | Ltr.(App) |
|-------------------|-----------------------------|-----------|
| OT-VO-136-MS/SS-A | 225 mm x 300 mm (9" x 12")  | 12        |
| OT-VO-136-MS/SS-B | 300 mm x 300 mm (12" x 12") | 22        |
| OT-VO-136-MS/SS-C | 300 mm x 450 mm (12" x 18") | 32        |
| OT-VO-136-MS/SS-D | 300 mm x 600 mm (12" x 24") | 43        |
| OT-VO-136-MS/SS-E | 375 mm x 500 mm (15" x 20") | 56        |
| OT-VO-136-MS/SS-F | 375 mm x 550 mm (15" x 22") | 60        |

**Optional:**

1. Digital Vacuum Gauge.
2. Microprocessor Based Touch Screen Display System.

**Note:** Without Vacuum Pump.

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: [info@opticstechnology.in](mailto:info@opticstechnology.in)

[opticstechnology@gmail.com](mailto:opticstechnology@gmail.com)

Website : [www.opticstechnology.in](http://www.opticstechnology.in)