



OPTICS TECHNOLOGY

SINCE 1978

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2013 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

VACUUM TRAY DRYER

RECTANGULAR CE



OT-VTD-138 MS



OT-VTD-138 SS

Vacuum Tray Dryer used in various industries under vacuum atmosphere to dry temperature-sensitive materials. It gives complete drying and vaporization of moisture of the product. Inside Chamber made of S.S. (304) & Outer Mild steel duly powder coated / S.S. Temperature controlled by Microprocessor based PID Controller, Temp. Stability $\pm 1^{\circ}\text{C}$, Temperature Range 50°C to 130°C . It Fitted with Vacuum Release Valve, with Toughened Glass Door, Provided with S.S. Tray. Tray size **16" x 32"**.

Model

No. of Tray

OT-VTD-138-MS/SS-A	4 No. S.S. Trays
OT-VTD-138-MS/SS-B	6 No. S.S. Trays
OT-VTD-138-MS/SS-C	8 No. S.S. Trays
OT-VTD-138-MS/SS-D	12 No. S.S. Trays
OT-VTD-138-MS/SS-E	24 No. S.S. Trays
OT-VTD-138-MS/SS-F	48 No. S.S. Trays

Optional:

1. Digital Vacuum Gauge
2. Microprocessor Based Touch Screen Display System with Data Logging
3. Vacuum Pump with Oil & Belt Motor System
4. Vacuum Gauge
5. Moisture Trap
6. Non-Return Valve

Note: Without Vacuum Pump.

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in

opticstechnology@gmail.com

Website : www.opticstechnology.in