



OPTICS TECHNOLOGY

SINCE 1978

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2018 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022



FREEZE DRYER (LYOPHILIZER)

LAB MODEL

CE



OT-FD-132

Condenser Temp. **(-) 40 °C / (-) 70 °C or better** (at ambient temp. of 20 °C) **PUF INSULATED, CFC Free Compressor**, Drying Manifold (S.S 304) ½" and ¾". Direct Drive Vacuum Pump 100Ltrs Minimum capacity, continuous display Temp. Indicator (LED Display), Digital Type Pirani Gauge, Housed in a Super Structure frame with wheel based.

Model No. Capacity

OT-FD-132-A 3 Ltrs With 6 port standard

OT-FD-132-B 5 Ltrs With 12 port standard

Optional

1. Digital Vacuum Indicator
2. Microprocessor based LCD display system for vacuum and temperature

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in

opticstechnology@gmail.com

Website : www.opticstechnology.in