



SINCE 1978

OPTICS TECHNOLOGY

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED
An ISO 10002:2018 CERTIFIED
An ISO 13485:2016 CERTIFIED

IEC 61010-1:2010
BS EN 12469:2000
ISO 29463-5:2022

PORTABLE AUTOCLAVE

CE

(Pressure Cooker Type)



OT-PAU-69 AL

ALUMINIUM BODY



OT-PAU-69/2 SS

Electrically heated, complete with stand single walled Aluminum Body/Stainless Steel 304 along with essential parts like pressure gauge, Safety valve etc.

Model No.	Size	Capacity (Approx.)
OT-PAU-69-AL/SS-11	Size 225 x 275 mm	11 Ltrs.
OT-PAU-69-AL/SS-21	Size 300 x 300 mm	21 Ltrs.
OT-PAU-69-AL/SS-25	Size 300 x 350 mm	25 Ltrs.
OT-PAU-69-AL/SS-28	Size 300 x 400 mm	28 Ltrs.
OT-PAU-69-AL/SS-38	Size 350 x 450 mm	38 Ltrs.
OT-PAU-69-AL/SS-39	Size 350 x 550 mm	39 Ltrs.

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in,

opticstechnology@gmail.com,

Website : www.opticstechnology.in