



OPTICS TECHNOLOGY



SINCE 1978

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED
An ISO 10002:2018 CERTIFIED
An ISO 13485:2016 CERTIFIED

IEC 61010-1:2010
BS EN 12469:2000
ISO 29463-5:2022

AUTOMATIC LEAK TEST APPARATUS

CE

(GMP MODEL)



OT-LTA-25/A

It is upgraded version of conventional leak test machine which adopts a user friendly and ergonomic method. With Digital Vacuum Gauge help in improving the accuracy and helps in releasing vacuum at a set time. With a microprocessor system which helps automatic process as per the sequence there by eliminating human/operator interventions. Vacuum level can be set table as per the requirement. It also has a provision to set hold time and retention time for better result.

Desiccator Diameter : 300mm (12" Dia)

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in,

opticstechnology@gmail.com,

Website : www.opticstechnology.in