



OPTICS TECHNOLOGY



SINCE 1978

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2018 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

AUTOMATIC BUTYRO REFRACTOMETER

(NO NEED OF CIRCULATING WATER BATH)



OT-ABR-16

Features:

- External Light Interference (ELI)
- Temperature Compensated scale at 40° C
- Automatic & Accurate (Accuracy: Butyro ± 0.5)
- Calibration with Standard liquid
- Large, easy to read display
- Display of both Butyro & RI scale
- **Automatic Temperature Compensated reading at 40°C**

| | |
|------------------------------|---|
| Range (Converted at 40°C) | Butyro 30.0 to 90.0 RI 1.4450 to 1.4850 |
| Minimum Indication | Butyro 0.1 RI 0.0001 |
| Accuracy | Butyro ± 0.5 (at 40°C) RI ± 0.0003 (at 40°C) |
| Temperature | 10°C-50°C Automatic Temp. Compensation |

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in,

opticstechnology@gmail.com,

Website : www.opticstechnology.in