



# OPTICS TECHNOLOGY



SINCE 1978

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010  
An ISO 10002:2018 CERTIFIED BS EN 12469:2000  
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

## AUTOMATIC MILK ANALYZER



OT-MA-24

Milk analyzer is designed especially for the Indian market. The device has a unique cleaning system with no need for Electricity and is very useful for work in field conditions. Correct Analysis of the components is carried out for 40 seconds. Stainless steel box, water-resistant Keypad, Protective Coating on the electronics prevent moisture problems.

### Features :

- No acid or other chemicals are used
- Ultrasonic Technology
- Very small quantity of milk
- Balloon pump with Manual cleaning function
- Self calibration without computer
- Low power consumption
- Large graphic display

### Parameters :

- FAT
- SNF
- Density
- Temp.
- Protein
- Lactose
- Added Water
- Salt

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: [info@opticstechnology.in](mailto:info@opticstechnology.in),

[opticstechnology@gmail.com](mailto:opticstechnology@gmail.com),

Website : [www.opticstechnology.in](http://www.opticstechnology.in)