



OPTICS TECHNOLOGY



SINCE 1978

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2018 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

AUTOMATIC REFRACTOMETER

(Digital Indication of Brix & Temperature)

Battery Operated



RANGE 0-80% BRUX

“Optics Technology” make Hand Refractometer (Brix Meter) is specially designed for measuring sugar content in an aqueous solution, soft drinks, sugar solutions, cutting liquid, fruit juice or milk and beverages.

Range	: 0 to 80% Brix
Resolution	: 0.1% Brix
Accuracy	: $\pm 0.2\%$ Brix
Temp. Resolution	: 0.1°C (0.1°F)
Temp. Compensation	: Automatic between 10°C and 40°C
Automatic Shut-Off	: After three minutes of inactivity

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in,

opticstechnology@gmail.com,

Website : www.opticstechnology.in