



SINCE 1978

# OPTICS TECHNOLOGY

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010  
An ISO 10002:2018 CERTIFIED BS EN 12469:2000  
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022



## DIGITAL DISSOLUTION RATE TEST APPARATUS WITH WATER BATH

(6 STAGE) CE



OT-DDRT-176/1

The unit is equipped with acrylic water bath controlled by digital temperature controller which is normally set at 37°C and duly fitted with pump which recirculates the water thus maintaining homogeneous temperature. The acrylic water bath is with digital timer can be set ranging from 1 to 999 minutes add it "BEEPS" at the end of the set time.

The unit is equipped with Digital setting of RPM as an optional which can be controlled from 25 to 150 RPM.

The equipment is duly fitted with motorized lifting to facilitate easy operation of the stirrers' assembly.

Digital dissolution rate test apparatus is equipped With

- Digital RPM Indicator - Range 25 to 150 RPM
- Digital Timer – Range 1 to 999 minutes
- Digital Temperature Indicator – Range 30°C to 50°C
- Body is available M.S. / S.S. 304

Available in :  
6 Stage &  
8 Stage

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: [info@opticstechnology.in](mailto:info@opticstechnology.in)

[opticstechnology@gmail.com](mailto:opticstechnology@gmail.com)

Website : [www.opticstechnology.in](http://www.opticstechnology.in)