



SINCE 1978

OPTICS TECHNOLOGY



Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2018 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

DIGITAL TEMPERATURE CUM HUMIDITY CONTROLLER

CE

Mains Operated



OT-DTEH-77

Features:

- 7 Segment LED Display
- High Accuracy
- Power Supply : 220 V AC $\pm$ 15% 50/60 Hz
- 1" Large LCD display wall mounted.
- Temperature Range : 0°C~50°C $\pm$ 0.5°C
- Humidity Accuracy : $\pm$ 3% RH
- Temperature Reading only in °C

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in,

opticstechnology@gmail.com,

Website : www.opticstechnology.in