



OPTICS TECHNOLOGY



SINCE 1978

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2018 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

DIGITAL THERMOHYGROMETER

CE

(Battery Operated)



OT-DTH-76

Features:

- Temperature Range : 0°C~50°C
- Temp. Accuracy : $\pm 1.0^{\circ}\text{C}$
- Humidity Range : 10%~99% RH
- Humidity Accuracy : $\pm 5\%$ RH, (40%~80%)
- Temp. Reading for $^{\circ}\text{F}$ & $^{\circ}\text{C}$
- MAX. and MIN. Temperature
- LCD Display Wall Mounted.

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in,

opticstechnology@gmail.com,

Website : www.opticstechnology.in