



OPTICS TECHNOLOGY



SINCE 1978

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2018 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

DIGITAL TURBIDITY / NEPHELO METER

CE



OT-DTM-184/2
DELUXE MODEL

	Normal Model OT-DTM-184/1	Deluxe Model OT-DTM-184/2
Range	0-1000 (in 2 Ranges) 0-200 & 0-1000	0 - 1000 NTU (in 4 Ranges) 0 - 2, 20, 200, 1000 NTU
Resolution	0.1 NTU in 0-200 1 NTU in 0-1000	0.01 NTU in 0-2 & 0-20 0.1 NTU in 0-200, 1 NTU in 0-1000

Accuracy	: ± 3% FS, ± 1Digit
Display	: 3½ Digit 7-Segment LED
Response Time	: Less than 10 seconds
Light Source	: 6V, 0.3Amp Tungsten Lamp
Detector	: Photodiode
Sample System	: 30 mm Clear Glass Test Tubes
Accessories	: Test tube Set of 5, Operation Manual

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in

opticstechnology@gmail.com

Website : www.opticstechnology.in