



OPTICS TECHNOLOGY

SINCE 1978

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2018 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

DIGITAL pH AND CONDUCTIVITY METER COMBINED (Table Top) CE



OT-DpHCM-108

	pH	CONDUCTIVITY	TEMPERATURE
Range	0-14	0 to 20mS/cm in 2 Ranges	0-100°C
Resolution	0.01	1μS/cm	0.1°C
Accuracy	± 0.01 pH ± 1 digit	± 0.5% ± 1 digit	± 0.2% ± 1 digit

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in

opticstechnology@gmail.com

Website : www.opticstechnology.in