



OPTICS TECHNOLOGY



SINCE 1978

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2018 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

DISINTEGRATION TEST APPARATUS

DOUBLE BASKET

CE



OT-DTA-173/2 SS

DOUBLE BASKET



DT WITH WATER BATH & DIGITAL TIME TOTALISER (DOUBLE BASKET)

DT WITH WATER BATH & DIGITAL TIME TOTALISER (DOUBLE BASKET)

The instrument is designed for accurate estimation of disintegration time of tablets/capsules. It is assembled in painted sheet metal cabinet with chemical resistant epoxy coating, fitted with superior quality constant speed sturdy motor. The guide discs and strokes are conforming to Pharmacopoeia Standards. (Body is made of M.S. Duly Powder Coated / SS 304 Body)

Specifications :

Displacement	: 55mm
RPM	: 30
Drive	: 1/8 HP
Heater	: 125 watt.
Power	: 230 V AC Single Phase

Optional:

- 1) Digital Time Totalizer
- 2) Water bath System Instead of Hot Plate with Thermostatics Temperature Controller.
- 3) Microprocessor based PID Controller (Only with Water bath model).

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in

opticstechnology@gmail.com

Website : www.opticstechnology.in