



SINCE 1978

OPTICS TECHNOLOGY

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2018 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

DISTILLATION APPARATUS

CE

(S.S. 304 WALL MOUNTING SYSTEM)

Deluxe Model Heavy Duty Argon Welded Model



OT-DA-WM-122

Model No.	Capacity
OT-DA-WM-122 A	4 Ltr/hr
OT-DA-WM-122 B	6 Ltr/hr
OT-DA-WM-122 C	8 Ltr/hr
OT-DA-WM-122 D	10 Ltr/hr

All parts of Distillation apparatus are made out of stainless steel. It consists of Main Chamber, Condenser Tube, Condenser Jacket, Vapour Baffle.

Optional:

1. Water level cutoff device (For Single Phase)
2. Water level cutoff device (For Three Phase)

Note:

We recommended to use water level cut-off device for safety of heaters & long life of instruments

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in

opticstechnology@gmail.com

Website : www.opticstechnology.in